

RDBLPR-520B ELISA WASHER



Characteristics:

- ☐ Large LCD display, intelligent software, easy to operate
- ☐ Strong washing function, bottom washing and two positions aspiration make completed washing effect
- ☐ 8-channel or 12 channels optional.
- ☐ Needle positions adjustable and up to 10 different needle positions pre-settable Pressure, volume & time adjustable, washing times, soap time pre-settable Pressure, volume & time adjustable, washing times, soap time pre-settable, whole plate or single strip reading
- ☐ Needle positions adjustable and up to 10 different needle positions pre-settable 5 positions of needles are adjustable, fit for multi-kind of plates.
- ☐ Self-calibration and trouble alarm function.
- ☐ Large memory to store up to 100 programs.



Specification:

Plate: 48/96 well and trips (U,V, Flat bottom)

Residual Volume: $\leq 2\mu\text{l}$

Dispensing volume: Adjustable from 50 to 900 μl

Wash cycles: 1-16 optional

Mode of wash: From 1 to 12, row and rank, selectable

Soak time: Up to 990 seconds adjustable`

Shake time: Up to 990 seconds adjustable.

Display: Large LCD

Channel: 3 channels for washing liquid, 1 channel for waste liquid

Memory: up to 100 programs can be stored

Interface: RS232

Dimension: 450(L) mm \times 310(W)mm \times 170 (H)mm

Net Weight: 7.5kg

Power Supply: AC220V \pm 10%, 50/60Hz; AC110V \pm 10%, 50/60Hz