

ICE FLAKING MACHINE



Features:

- Fully Automatic microprocessor control
- Compact Design and low maintenance
- Stainless Steel construction
- Continuous Ice Flakes Output
- CFC free compressor
- Noiseless operation
- Overload protection
- Low water level detection
- Wheels for easy mobility

Specifications:

MODEL	RDBL-IF-30	RDBL-IF-50	RDBL-IF-85	RDBL-IF-100
Capacity Of Ice Flake(Kg Pr24 Hrs)	30kg	50kg	85kg	100kg
Storage Bin	Abs Plastic			
Storage Capacity	10kg	15kg	25kg	
Cooling Method	Air			
Working Temperature	Room temp should be below 30°C			
Construction	Stainless steel			
Noise Level	≤52dB			
Water Supply	¾ inch BSP			
Refrigerant	R134A			
Input Power	180w	280w	420w	
Dimension(Mm)	30000x493x547	380x543x722	549x611x833	
Power Supply	220v /50HZ			
Gross Weight	36kg	45kg	62kg	64kg

Applications:

Research Institute, Hospitals, Laboratories, Biotechnology & Universities.