

## RDBL-MINI-6000 MINI FUGE



### Description:

The mini centrifuge series MINI-6K is perfect designed with novel and unique appearance. It is pretty and skillful, equipped with 2 types of rotors and several tube carriers for 2.0ml/ 1.5ml / 0.5ml / 0.2ml tubes and 0.2ml x 8 PCR tube rank. It is designed for feeling. Open/close the cover will automatically start/stop the centrifuge. Built-in timing function. Humanism is fully displayed in its transparent roundness cover, several rotors, and tube carriers for choice.

## Technical Parameters:

Model	Mini-4K	Mini-6K	Mini-7K	Mini-10K	Mini-10K+
Speed	4000rpm	6000rpm	7000rpm	10000rpm	3000rpm - 10000rpm (step 1000rpm)
Centrifugal Force	1200g	3000g	3400g	7500g	1000g - 7500 g
Sample Capacity	8 x2.0ml/1.5ml/0.5ml/0.2ml tube (additional 0.5ml and 0.2ml tube carriers) 2 x 8 x 0.2ml PCR tube strip		8 x2.0ml/1.5ml/0.5ml/0.2ml tube (additional 0.5ml and 0.2ml tube carriers )		8 x 2.0ml/1.5ml/0.5ml/0.2ml tube (additional 0.5ml and 0.2ml tube carriers)2 x 8 x 0.2ml PCR tube strip(tube strip only used when speed≤6000rpm)
Time Range	1s - 99min59s				1s - 999min
Noise	≤45 dB		≤55dB		
Dimension	176x156x121 (mm)				
Net Weight	≤1.5kg				
Power	220V/110V 50-60HZ				