

MINI VERTICAL ELECTROPHORESIS 2 GEL



Technical Specifications:

Specifications	RDPRO-2 Mini Vertical Blot (Gel Quantity: 2 gels)
Gel quantity	1-2PCS
Gel thickness	1.0mm (standard)
Combs	10 wells/15 wells(standard)
Maximum volume of sample	10 wells:44μl 15 wells:26μl
Buffer Capacity	700ml
Running time	It is related to the power setting, which should not be set too large. Otherwise the strip(band) will be discrete.
Running gel sizes	101*73mm
Removable UV-Transparent Gel Tray	-
Maximum voltage limit	-
Maximum power limit	-
Maximum buffer temperature	-
Working conditions--temperature	-
Working conditions--humidity	-
Dimension (L×W×H)	180×130×160mm
Net Weight	2kg