

RDBL-PH-10 PH METER



Specification of pH meter:

- pH Range -2.00 to 16.00Ph
- pH Resolution ± 0.01 pH.
- pH Accuracy ± 0.01 pH.
- mV Range 0 to +1999mV.
- mV Accuracy ± 0.01 mV.
- Temperature Range 0 to 100o C.



- Temperature Accuracy $\pm 1^{\circ}\text{C}$.
- Calibration Points 3 Points.
- pH Buffer Options NIST (pH 4.00/6.86/9.18) or (pH4.01/ 7.00/10.00)
- Temperature Compensation 0 to 100°C , Manual or Automatic.
- Power Requirements DC 9V, using AC Adaptor 220 VAC, 50 Hz.
- Dimensions 210 (L) mmX205 (W) mm X 65 (H) mm.
- Weight 1.75kg