

RDBLWB-203PWEIGHING BALANCE



Description:

The visual appearance of the product is blunt into the simple concept, stylish and Generous. Specially makes the fine pouring cast- aluminium outer covering, so as to enhance the balance of the anti-static anti-jamming capability.

The panel with five operational buttons, Sliding glass windscreen provide big space, insightful vision and operate easily. We use high precision magnet sensor to keep the weighing faster and stable.

Magnet sensor keep weighing stable and fast

Full range tare, Counting, Units choose

RS 232 Interface can connect the computer and printer



TECHNICAL PARAMETERS:-

Model	RDBLWB-203P
CAPACITY	200gm
CLASS	1
RESOLUTION	0.001g (1mg)
REPEATABILITY	0.001g
LINER 0.002g	
STABLE TIME	≤3s
OP.TEMP	10°C~40°C
PANSIZE	Φ80mm
DEMENSION	415 x 405 x 295 (L×W×H) mm
POWER	Power Line
CAL.METHOD	EXTERNAL