

## RDBL- MIX1000 VORTEX MIXER



### Description:

MIX1000 Mini Vortex Mixer is designed for “touch” mixing a variety of test tubes. With a 4.5mm orbit and adjustable speed of 1000-2800rpm, the mini vortex mixer instantly vortexes even the largest samples, including nearly full 50ml tubes. The small footprint allows it to be used on even the most crowded bench

### Features:

Shaking speed: 1000-3500rpm.

Powerful and reliable motor designed for daily laboratory use. “Touch” pressure activated operation.

Accept a variety of test tubes, including the tubes up to 50ml. Mini-sized, allowing it to fit on even the most crowded bench



## Parameters:

Project	Parameters:
Speed	1000-3500rpm
Type of movement	Orbital
Max test tube size	Φ30mm
Voltage	AC100-240V, 50/60Hz
Dimension(mm)	125×115×76
Type of motor activation	Touch pressure activated operation
Orbit	4.5mm
Power	24W
Weight(kg)	0.55