

RDBLWJJ2204- WEIGHING BALANCE



Description:

The visual appearance of the product is blunt into the simple concept, stylish and Generous. Specially makes the fine pouring cast- aluminium outer covering, so as to enhance the balance of the anti-static anti-jamming capability. The panel with Five operational buttons, sliding glass windscreen provide big space, insightful vision And operate easily. We use high precision magnet sensor to keep the weighing moreFast and stable.

Magnet sensor keep weighing stable and fast

Full range tare, Counting, Units choose

RS 232 Interface can connect the computer and printer



TECHNICAL PARAMETERS:-

| | |
|---------------|------------------------|
| Model | RDBLWJJ2204 |
| CAPACITY | 220gm |
| CLASS1 | |
| RESOLUTION | 0.0001g (0.1mg) |
| REPEATABILITY | ±0.0002g |
| LINER | ±0.0003g |
| STABLE TIME | ≤3s |
| OP.TEMP | 17.5°C~22.5°C |
| PANSIZE | Φ80mm |
| DEMENSION | 340*215*350 (L×W×H) mm |
| POWER | Power Line |
| CAL.METHOD | EXTERNAL |
| INTERFACE | RS232 |