

ULTRASONIC BATH/CLEANER



Features:

- 0~30 mins digital timer;
- 20~80°C temperature digital control;
- 2.0mm thickness stainless steel tank;
- All models with drainage system and heat cooling fan;
- 40~100% ultrasonic power adjustable;
- Strong mother board and allow for long time working,
- Welding tank.



TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer wood packaging size	Total G.W (KG)
RDBPL-060AL	14	360	500	300x240x200	400x300x380	430x380x450	14.0	Carton packaging	14.0
RDBPL-080AL	20	480	500	330x300x200	430x350x380	460x430x450	18.0	Carton packaging	18.0
RDBPL-100AL	30	600	500	500x300x200	600x350x380	670x485x470	26.0	670x485x470	33.0
RDBPL-120AL	38	720	1000	500x300x250	600x350x480	710x460x630	31.0	710x460x630	39.0
RDBPL-180AL	58	1080	1000	550x350x300	650x400x580	710x450x720	42.0	710x450x720	51.0
RDBPL-240AL	78	1200	1000	550x400x350	650x450x620	720x540x830	51.0	720x540x830	68.0
RDBPL-360AL	130	2160	1000	650x450x450	750x500x720	810x560x910	75.0	810x560x910	88.0

ULTRASONIC BATH/CLEANER



Features:

- Double frequency 28KHZ or 40KHZ can be selected for working;
- Degas function;
- 0~99 mins digital timer;
- 20~80°C temperature digital control;
- 6L upper machine installed drainage system;
- 10L upper machine installed heat dissipation fan.



TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-10AD	2	80	100	150x140x100	190x170x220	250x250x320	3.0	540x540x650(8pcs)	25.0
RDBPL-20AD	3	120	100	240x135x100	270x170x240	340x250x320	4.0	790x380x690(6pcs)	25.0
RDBPL-30AD	6	180	300	300x150x150	330x180x310	430x280x420	6.0	860x460x820(6pcs)	38.0
RDBPL-40AD	10	240	300	300x240x150	330x270x310	430x370x430	8.0	760x460x780(4pcs)	34.0
RDBPL-60AD	15	360	400	330x300x150	360x330x310	460x410x430	10.0	880x490x790(4pcs)	42.0
RDBPL-80AD	22	480	600	500x300x150	550x330x310	660x440x430	13.0	670x460x790(2pcs)	28.0
RDBPL-100AD	30	600	800	500x300x200	550x330x360	660x440x480	15.0	670x460x890(2pcs)	32.0

ULTRASONIC BATH/CLEANER



Features:

- 0~30 mins mechanical timer;
- 20~80°C temperature mechanical control;
- 6L upper machine installed drainage system;
- 10L upper machine installed cooling fan.

TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-08	1.3	60	50	150x140x65	180x170x140	260x260x260	2.5	540x540x600(8pcs)	20.0
RDBPL-10	2	80	50	150x140x100	180x170x180	260x260x290	3.0	540x540x650(8pcs)	25.0
RDBPL-20	3	120	100	240x135x100	270x170x210	350x260x330	4.0	790x380x690(6pcs)	25.0
RDBPL-30	6	180	200	300x150x150	330x180x270	430x280x380	6.0	860x460x820(6pcs)	38.0
RDBPL-40	10	240	250	300x240x150	330x270x270	430x370x380	8.0	760x460x780(4pcs)	34.0
RDBPL-50	14	300	400	300x240x200	330x270x320	430x370x430	9.0	760x460x880(4pcs)	42.0
RDBPL-60	15	360	400	330x300x150	360x330x270	460x430x390	10.0	880x490x790(4pcs)	42.0
RDBPL-70	19	420	600	330x300x200	360x330x320	460x430x430	11.0	880x490x890(4pcs)	48.0
RDBPL-80	22	480	600	500x300x150	550x330x270	630x430x390	13.0	670x460x790(2pcs)	28.0
RDBPL-100	30	600	800	500x300x200	550x330x320	630x430x430	15.0	670x460x890(2pcs)	32.0

ULTRASONIC BATH/CLEANER



Features:

- 0~30 mins digital timer;
- 20~80°C temperature digital control;
- 6L upper machine installed drainage system;
- 10L upper machine installed heat cooling fan.

TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-08A	1.3	60	50	150x140x65	180x170x140	260x260x260	2.5	540x540x600(8pcs)	20.0
RDBPL-10A	2	80	50	150x140x100	180x170x180	260x260x290	3.0	540x540x650(8pcs)	25.0
RDBPL-20A	3	120	100	240x135x100	270x170x210	350x260x330	4.0	790x380x690(6pcs)	25.0
RDBPL-30A	6	180	200	300x150x150	330x180x270	430x280x380	6.0	860x460x820(6pcs)	38.0
RDBPL-40A	10	240	250	300x240x150	330x270x270	430x370x380	8.0	760x460x780(4pcs)	34.0
RDBPL-50A	14	300	400	300x240x200	330x270x320	430x370x430	9.0	760x460x880(4pcs)	42.0
RDBPL-60A	15	360	400	330x300x150	360x330x270	460x430x390	10.0	880x490x790(4pcs)	42.0
RDBPL-70A	19	420	600	330x300x200	360x330x320	460x430x430	11.0	880x490x890(4pcs)	48.0
RDBPL-80A	22	480	600	500x300x150	550x330x270	630x430x390	13.0	670x460x790(2pcs)	28.0
RDBPL-100A	30	600	800	500x300x200	550x330x320	630x430x430	15.0	670x460x890(2pcs)	32.0

ULTRASONIC BATH/CLEANER



Features:

- 0~30 mins digital timer;
- 20~80°C temperature digital control;
- 40~100% ultrasonic power adjustable;
- 6L upper machine installed drainage system;
- 10L upper machine installed cooling fan.



TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-20AL	3	50~120	100	240x135x100	270x170x210	350x260x330	4.0	790x380x690(6pcs)	25.0
RDBPL-30AL	6	70~180	200	300x150x150	330x180x270	430x280x380	6.0	860x460x820(6pcs)	38.0
RDBPL-40AL	10	100~240	250	300x240x150	330x270x270	430x370x380	8.0	760x460x780(4pcs)	34.0
RDBPL-50AL	14	120~300	400	300x240x200	330x270x320	430x370x430	9.0	760x460x880(4pcs)	42.0
RDBPL-60AL	15	150~360	400	330x300x150	360x330x270	460x430x390	10.0	880x490x790(4pcs)	42.0
RDBPL-70AL	19	170~420	600	330x300x200	360x330x320	460x430x430	11.0	880x490x890(4pcs)	48.0
RDBPL-80AL	22	200~480	600	500x300x150	550x330x270	630x430x390	13.0	670x460x790(2pcs)	28.0
RDBPL-100AL	30	240~600	800	500x300x200	550x330x320	630x430x430	15.0	670x460x890(2pcs)	32.0

ULTRASONIC BATH/CLEANER



Features:

- 0~99 mins digital timer;
- 20~80°C temperature digital control;
- 40~100% ultrasonic power adjustable(digital control power adjustabel);
- 6L upper machine installed drainage system;
- 10L upper machine installed cooling fan.



TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power(W)	Heater power(W)	Tank size (LxWxH)mm	Machine size (LxWxH)mm	Unit Packaging size(LxWxH) mm	G.W. (KG)
RDBPL-1340AL	4	180	100	300x155x100	330x180x240	410x250x340	5.0
RDBPL-1360AL	6	180	300	300x155x150	330x180x310	440x270x400	6.5
RDBPL-1410AL	10	240	300	300x240x150	330x270x310	430x350x410	9.5
RDBPL-1615AL	15	360	400	330x300x150	360x330x310	460x400x410	11.0
RDBPL-1820AL	22	480	500	500x300x150	550x330x310	670x440x420	14.5
RDBPL-11030AL	30	600	500	500x300x200	550x330x360	670x440x470	16.5

ULTRASONIC BATH/CLEANER



Features:

- 0~30 mins mechanical timer;
- 6L upper machine installed drainage system;
- 10L upper machine installed cooling fan.

TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-08T	1.3	60	None	150x140x65	180x170x140	260x260x260	2.5	540x540x600(8pcs)	20.0
RDBPL-10T	2	80	None	150x140x100	180x170x180	260x260x290	3.0	540x540x650(8pcs)	25.0
RDBPL-20T	3	120	None	240x135x100	270x170x210	350x260x330	4.0	790x380x690(6pcs)	25.0
RDBPL-30T	6	180	None	300x150x150	330x180x270	430x280x380	6.0	860x460x820(6pcs)	38.0
RDBPL-40T	10	240	None	300x240x150	330x270x270	430x370x380	8.0	760x460x780(4pcs)	34.0
RDBPL-50T	14	300	None	300x240x200	330x270x320	430x370x430	9.0	760x460x880(4pcs)	42.0
RDBPL-60T	15	360	None	330x300x150	360x330x270	460x430x390	10.0	880x490x790(4pcs)	42.0
RDBPL-70T	19	420	None	330x300x200	360x330x320	460x430x430	11.0	880x490x890(4pcs)	48.0

ULTRASONIC BATH/CLEANER



Features:

- LED display;
- Big volume: 2500ML;
- Five interval time: 90s/180s/280s/380s/480s;
- Heating function
- Available for cleaning Jewelry, watches, rings, necklace, earrings, diamond, glasses, sunglasses, optical lenses, dentures, combs, tooth brushes, shaver head.



TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-4820	2.5	60	100	250x150x80	290x223x185	345x250x230	3.2	520x350x480(4pcs)	14.0

ULTRASONIC BATH/CLEANER



Features:

- Ultrasonic power: double power can be choosed 30W / 50W;
- Timer adjustable: 1-60 min (30W), 1-30min (50W);
- Working frequency: 40KHZ

TECHNICAL SPECIFICATIONS:-

Model	Capacity (L)	Ultrasonic power (W)	Heater power (W)	Tank size (LxWxH) mm	Machine size (LxWxH) mm	Unit Packaging size (LxWxH) mm	G.W. (KG)	Outer Carton size (LxWxH) mm	Total G.W (KG)
RDBPL-968	0.6	50	None	160x80x38	185x145x100	215x170x130	0.72	440x360x530 (16pcs)	12.5